

Sputocol™ Sputum Collection System



Sputocol™ Sputum Collection System

Our patented Sputocol™ is a self-contained device for collecting, transporting and processing sputum. The collection tube is not in contact with the patient or the laboratory personnel during the collection process until it is capped.

Patients expectorate into the large funnel and reclose the lid. The inside edge of the funnel brings the specimen directly into the collection tube and prevents contact with the tube rim.

The bottom lid detaches and functions as a capping device so that personnel can then prepare the specimen for transport without contacting a contaminated surface area.

Sputocol™ is available with a conical bottom tube or a freestanding base tube. A patient label is included with every unit for convenience.

Material: Polypropylene collection tube, Polyethylene cap.

For more information, to request a sample or to contact your sales representative, please contact us:

800.421.6261
info@evergreensci.com
www.evergreensci.com

Catalog Number	Description	Packaging
221-3995-G80	Sputocol System, Freestanding 50mL Tube	Indv. Pack, 72/case
221-3996-G80	Sputocol System, Conical 50mL Tube	10/bag, 60/case

